



CryoCube

- A New First in Market low cost dry shipper, which combines long hold times with competitive pricing without sacrificing sample viability in transportation.
- Features include – Safe shipping in any orientation, suitable for international transportation, weighs less than 4.1Kg (9lbs) when charged.



Doble Qwick Series

- All of the quality, reliability and robust build features of the standard Doble series but these units charge in less than two hours, allowing for same-day vapour shipping.
- Features include – Hold times of up to 20 days (vapour). Low liquid nitrogen consumption. Dual function, as once delivered to the point of use, units can be filled with nitrogen and utilised for long term sample storage.



HE Series Freezers

- Highest specification freezer available in the marketplace, with -190°C certified temperature below the lid.
- Features include – 94,000 vial capacity, state of the art electronic controller unit, hold time (from last timed fill) up to 20 days thus ensuring sample viability in cases of either a catastrophic failure in the nitrogen supply system or other unforeseen emergencies. Lowest like for like running costs against any competitor unit.



Research Dewars

- Designed for general liquid nitrogen transportation and snap freezing of samples.
- Features include – Capacities from 0.5 to 6 litre, integrated carry handle, high strength stainless steel construction, five year vacuum warranty.



Safety Training Courses

- Full range of gas safety training courses, available either on-line (£49.50) or workshop based (P.O.A). Subjects include: Using Liquid Nitrogen Safely, Cryogenic Gases Safety Awareness Workshop, Laboratory Gases Safety Awareness Workshop and Liquid Nitrogen Decanting Practical.

For more information on any of the above products go to www.cryoservice.co.uk or contact the Sales Office on **01905 758200**.